


# Patient Dose in DXA

Noirin Sheahan, PhD  
Medical Physics and Bioengineering,  
St. James's Hospital


## Outline of Presentation

1. Relevant anatomical considerations
2. Dose and Risk in DXA  
Adults  
Paediatrics
3. Foetal Dose during Pregnancy
4. Summary


### 1. Relevant Anatomical Considerations



### Irradiated Area




Radiation scattered outside of beam area contributes to dose:




### Radiosensitive Tissues in Irradiated Area

#### **Weighting Factors**

Bone marrow 0.12




### Radiosensitive Tissues in Irradiated Area



#### **Weighting Factors**

Bone marrow	0.12
Stomach	0.12


### Radiosensitive Tissues in Irradiated Area



#### **Weighting Factors**

Bone marrow	0.12
Stomach	0.12
Colon	0.12


### Radiosensitive Tissues in Irradiated Area



#### **Weighting Factors**

Bone marrow	0.12
Stomach	0.12
Colon	0.12
Liver	0.05


### Radiosensitive Tissues in Irradiated Area



#### **Weighting Factors**

Bone marrow	0.12
Stomach	0.12
Colon	0.12
Liver	0.05
Bladder	0.05


### Radiosensitive Tissues in Irradiated Area



#### **Weighting Factors**

Bone marrow	0.12
Stomach	0.12
Colon	0.12
Liver	0.05
Bladder	0.05
Gonads	0.20
	(.05)


### Radiosensitive Tissues in Irradiated Area




#### **Weighting Factors**

Bone marrow	0.12
Stomach	0.12
Colon	0.12
Liver	0.05
Bladder	0.05
Gonads	0.20
	(.05)
Skin	0.01
Bone Surface	0.01

**Variation of Dose with Depth –Under Couch Tube**  
e.g. Lunar DPX/Prodigy & Hologic




**Variation of Dose with Depth – Over Couch Tube**  
e.g. Lexxos & Lunar Expert





**2. Reported Dose Estimates in DXA**






Estimation of Effective Dose:



**Effective Dose (uSv) – Adults**  
(Spine plus dual hip – Standard Patient)

	Pencil	Fan	Cone
<b>Lunar</b> (DPX/Prodigy, std patient mode)	0.5	2	
<b>Hologic</b> (QDR1000/4500, Array mode)	6	32	
<b>Lexxos</b> (small/std patient)			18
<i>Lunar Expert:</i>		170	
<i>For comparison:</i>	<i>Dental X-Ray 1 – 8uSv</i> <i>Chest XRay 30-50uSv</i>		





**Effective Dose - Child compared to Adult**

Hologic Fan Beam (QDR 4500A – Fast Array Mode)			
	10 years		5 years
	Scaled		Scaled
Spine	x1		x2
Hip (male)	x2		x3
Hip (female)	x1		x1

*Thomas et al 2005*



Hologic Fan Beam (QDR 4500A – same mode for adult & child)				
	10 years		5 years	
	Scaled Default		Scaled Default	
Spine	x2	x3	x2	x4
Hip	x2	x2	x2	x3

*Blake et al 2006*


### Effective Dose - Child compared to Adult

Lunar Pencil Beam	10 years	5 years
	Scaled	Scaled
Spine	0.9	1.3


*Njeh et al 1997*



### 3. Foetal Dose during Pregnancy



Foetus receives scatter dose even in early pregnancy



### Foetal Dose from Hologic QDR 1000 Pencil Beam Scanner (Under-couch system)

	PA Spine (uSv)	Femur (uSv)	Spine+Dual-femur (uSv)
First trimester	2	3	7
Second trimester	3	1	6
Third trimester	5	1	7

*Damilakus et al 2002*


Foetal Natural Background Dose about 3uSv / day  
Foetal dose from AP Pelvis radiograph 800-1600uSv

No published reports but expect:

- similar figures or less for Lunar DPX / Prodigy
- about 5 times higher for Hologic Fan Beam


*Foetal dose from over-couch tube systems not yet reported ...*

*... expect higher dose to foetus*




*Under-couch systems – foetus protected by Sacrum*

*n.b. Damilakis et al assumed foetus outside main beam during first trimester. But ... narrow-frame woman?*




**Scale field wherever possible to exclude uterus**

#### Foetal Dose variation with Depth in Abdomen



#### Foetal Dose variation with Depth in Abdomen



#### Patient Dose in DXA - Summary

- Factor of 60 variation in dose (0.5 – 30uSv) depending on technology, supplier, model, mode (excluding Lunar Expert – high dose model)
- Similar to Dental - Chest X-Ray dose range.
- Few days natural background at most.
- Dose in Children may be up to 3 times adult dose – depending on technology, supplier, mode etc
- Foetal dose about 7uSv (pencil beam) – few days natural background at most.

#### Keeping Patient Dose Low

- Many radiosensitive organs in or near scan area, so
  - Choose model which allows user to vary scan area
  - Vary scan area according to patient size
- For paediatric work, choose model which incorporates low-dose paediatric exposure mode
- Empty / fill bladder as appropriate before DXA in early pregnancy or possible pregnancy
- Scan femur only (not spine) in late pregnancy?